HYCISUN[®]

Anbio Biotech

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REF: A6061213

Carton Size: 54 X 35 X 36 cm
G.W.: KG
N.W.: KG
Qty: 300 tests / CTN

Anbio(Xiamen)Biotechnology Co.,Ltd.
No.2016,Wengjiao West Road,
Xinyang Street,Haicang
District,361026 Xiamen,Fujian,China.

Sunbeam International GmbH
Schumanstr.12, 52146, Würselen,
Germany

Lotus Nl B.V.
Koningin Julianaplein 10,1e
Verd,2595AA,The Hague,Netherlands.

4 260676 530614

HYGISUN[®]

Anbio Biotech

**Rapid COVID-19 Antigen Test
(Colloidal Gold)**

REF: A6061213

Carton Size: 54 X 35 X 36 cm
G.W.: KG
N.W.: KG
Qty: 300 tests / CTN

Anbio(Xiamen)Biotechnology Co.,Ltd.
No.2016,Wengjiao West Road,
Xinyang Street,Haicang
District,361026 Xiamen,Fujian,China.

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Rapid COVID-19 Antigen Test (Colloidal Gold)

Saliva test for personal use

Packungsinhalt: **1** TEST
Package: **1** TEST

Anbio (Xiamen) Biotechnology Co., Ltd.
No.2016, Wengqiao West Road, Xinyang
Street, Huzhang District, Xiamen, Fujian,
China

CE REP

Leibes NL B.V.
Koningin Juliangelaan 10, 1e Vent, 2385AA,
The Hague, Netherlands

Sunbeam International GmbH
Schumannstr. 12, 52146, Würselen, Germany

Packungsinhalt: **1** TEST
Package: **1** TEST

Anbio (Xiamen) Biotechnology Co., Ltd.
No. 2016, Wengjiao West Road, Xinyang
Street, Haicang District, Xiamen, Fujian,
China

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

**Rapid COVID-19 Antigen Test
(Colloidal Gold)
Saliva test for personal use**


Packungsinhalt: **1 TEST**

Lot No: 20210101
Exp. Date: 20220101


Anbio (Xiamen) Biotechnology Co., Ltd.
Address: No. 100, Xiang'an Road, Xiamen, Fujian, China
Tel: +86 592 5511111
Fax: +86 592 5511112
Email: info@anbio.com
Website: www.anbio.com



Authorization Letter

Subject : Letter of Authorization

We, Anbio (Xiamen) Biotechnology Co.,Ltd. with detail address at No.2016,Wengjiao West Road,Xinyang Street, Haicang District, Xiamen,Fujian,P.R.China 361026, manufacturer hereby authorize Sunbeam International GmbH, Schumanstr.12, 521 46 Wurselen, Germany, as our authorized distributor to distribute, market our below products in Germany.

NO.	Product name	CE number
1	Rapid COVID-19 Antigen Test (Colloidal Gold)	NL-CA002-2020-53206
2	Neutralizing Antibody of SARS-CoV-2 Test (Immuochromatographic Assay)	NL-CA002-2020-54641
3	Rapid Neutralizing Antibody of SARS-CoV-2 Test (Colloidal Gold)	NL-CA002-2020-54641
4	Immuofluorescence Analyzer	NL-CA002-2020-50112

This letter of Authorization is valid from the date of issue for a period of two years.

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Feb 05, 2021

To whom it may concern,

We, Anbio(Xiamen) Biotechnology Co.,Ltd recently has launched our new version COVID19 Antigen test by nasal swab, throat swab, saliva 3 specimens which could greatly improve the specificity of the test by saliva specimen.

In addition, we hereby declare that our COVID19 Antigen test is effective to mutant strain in British and South Africa. Please see below details.

1. The British and South Africa mutant, commonly known as the B.1.1.7 subtype, has been named VOC-202012/01 novel coronavirus. The most important mutations in this strain relative to the common strain are three mutations in the RBD region K417N, E484K, and N501Y. These three mutations make it easier for the virus to bind to the ACE2 receptor in the human body, making it easier for the virus to spread.

2. Why our product is still effective: a. What we detected is mainly the N protein of novel coronavirus, which is relatively conservative, and this British & South Africa mutant strain is mainly the gene encoding RBD protein; b. The mutant region of the British & South Africa mutant strain is not within our detection range, so our detection reagent was still effective in detecting the antigen of the mutant strain, but it also means that our reagent could not distinguish the mutant strain from the common strain;

3. Simulated data:

3.1. Detection of recombinant full-length RBD protein of mutant strains:

Detection of recombinant full-length RBD protein of mutant strains of K417N,E484K,N501Y Mutated RBD Protein

mutated recombinant protein	common strain	Sensitivity of national reference
800pg/mL:G4	800pg/mL:G4	S1:G4 (6800TCID50)
400pg/mL:G3	400pg/mL:G3	S2:G4 (3400TCID50)
200pg/mL:G2	200pg/mL:G2	S3:G3 (1700TCID50)
100pg/mL:G1	100pg/mL:G1	S4:G2 (850TCID50)
50pg/mL:G0	50pg/mL:G0	S5:G0 (425TCID50)

Note: color darker than G0 can be judged positive

3.2. Detection of animal model samples: detection of S protein in Clara cells of respiratory tract infected in mouse model:


concentration (pg/mL)	1600	800	400	200	100	50
results	G5	G4	G3	G2	G0	G0

Your sincerely

JieLi Zhang

R&D Director

Anbio Biotechnoloty

DECLARATION OF CONFORMITY

According Directive 98/79/EC on In Vitro Diagnostic Medical Devices, Annex III.

Manufacturer: Anbio (Xiamen) Biotechnology Co.,Ltd.

Address: No.2016, Wengjiao West Road, Xinyang Street, Haicang District,361026 Xiamen, Fujian, China.

European Representative: Lotus NL B.V.

Contact person: Peter **E-mail:** peter@lotusnl.com

Address: Koningin Julianaplein 10,1e Verd, 2595AA,The Hague, Netherlands.

In Vitro Diagnostic Directive:

- **Rapid COVID-19 Antigen Test (Colloidal Gold)**

Category: Others.

Conformity assessment route: Declaration of Conformity IVDD Annex III

Applicable Standards:

ISO 13485:2016

ISO 14971:2019

EN ISO 18113-1:2011

EN ISO 18113-2:2011

EN ISO 18113-3:2011

EN 13641:2002

ISO 15223-1:2016

EN 13612:2002

ISO 23640:2015

EN 62366-1:2015

We, the manufacturer, herewith declare with sole responsibility that our product/s mentioned above meet/s the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In Vitro Diagnostic Medical Devices.

We agree to develop, implement and maintain a documented post-production monitoring process.

Signed on:04/09/2020

Place:Xiamen,China

Name of authorized signatory: *Daming Wang*

Position held in the company: General Manager

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
A handwritten signature in blue ink, consisting of a series of loops and curves, positioned above the name Dr. M.J. van de Velde.

Dr. M.J. van de Velde

BUREAU VERITAS
Certification

Anbio (Xiamen) Biotechnology Co., Ltd

No.2016, Wengjiao West Road, Xinyang Street, Haicang District, Xiamen City,
P.R. China

Certified site:

NO.2016, WENGJIAO WEST ROAD, XINYANG STREET, HAICANG DISTRICT, XIAMEN CITY,
P.R. CHINA

*Bureau Veritas Italia S.p.A. certifies that the Management System of the
above organisation has been audited and found to be in accordance
with the requirements of the management system standards detailed below*

EN ISO 13485:2016

Scope of certification

Design and manufacture of fluorescence-based immunoassay reagent and devices, as an aid in clinical assessment of cardiovascular, gastric, inflammation, diabetic and infectious diseases detection, as well as hormone, vitamin testing. Design and manufacture of IVD reagents as an aid in clinical assessment of blood type testing.

Certificate awarded in conformity with the requirements of ACCREDIA DT 02-DC Rev.00

Original cycle start date: **22/06/2020**

Expiry date of previous cycle: **n.a.**


Certification / Recertification Audit date: **17/05/2020**

Certification / Recertification cycle start date: **22/06/2020**

Subject to the continued satisfactory operation of the organization's Management System, this certificate expires on: **21/06/2023**

Certificate No. - Version: **IT298645-1**

Revision date: **22/06/2020**

GIORGIO LANZAFAME - Local Technical Manager

Certification body address:
Bureau Veritas Italia S.p.A., Viale Monza, 347 - 20126 Milano, Italia

SGQ N° 009A

Member of the Accreditations of Mutual Recognition EA, IAF and ILAC
Signatory of EA, IAF and ILAC Mutual Recognition Agreements

Further clarifications regarding the scope of this certificate and the applicability of the management system requirements may be obtained by consulting the organisation.
To check this certificate validity please refer to the website www.bureauveritas.it

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